


Orthodontic Auxiliaries

*Smiles made daily
in New England, USA*

*Contract Manufacturing with
Excellent R&D Assistance*



established 1947



Springs

36" long stainless steel open wound and closed wound springs

Stainless steel rotation, upright, "T" springs and anterior torquing springs

Nickel titanium open or closed coil in 7", 15" or 21" lengths

15mm long nickel titanium open coil spring with closed ends in light, medium and heavy force

Nickel titanium closed coil extension springs with eyelets in light, medium and heavy force in 9mm, 12mm, 14mm, 16mm and 18mm lengths

Micro Dental White coated nickel titanium springs available in open wound, closed wound and extension springs with eyelets

Nickel titanium stop wound springs (open and closed molar distalizing spring) in 7" and 15" lengths

7" lengths packaged in tubes, longer lengths on bobbins and shorter springs packaged in zip-lock bags

Ligatures

Made from high quality annealed stainless steel

Long standard preformed ligatures available with or without PTFE coating


Short, twisted ligatures available with or without PTFE coating

Coating colors available in almond, white, blue, green, pink, and purple

Flat short twisted Kobayashi ligatures available with almond coating or without coating

Long Kobayashi ligatures uncoated

Packaged in tubes or zip-lock bags





316LVM Flattened Multi-Strand

Soft wire made from high quality braided 316LVM stainless steel (patent pending)

Design offers improved oral hygiene; chain designs enable food trap, whereas flattened multi-strand has no large gaps

316LVM offers a lower modulus of elasticity vs. 304 stainless steel

Low carbon 316 for improved corrosion resistance. Vacuum melted for fewer impurities making it even more resistant to tarnish, discoloration, and other forms of damage

316LVM stainless steel is commonly used for implantable medical devices

Improved formability with low profile malleable 316 allows for improved bondable surface vs. braided/chain or solid bar designs offering optimal flex level for chewing

Wire will not fracture or fray

6" length packaged 10 per tube

60" length packaged on bobbin

Wireforms

Omega adjustment stops, packaged 100 per tube

Brass separating rings, .020 and .025, packaged 100 per zip-lock bag

Steiner wedges, small and large, packaged 10 per tube

Palatal bars, distal loop, or mesial loop from 35mm to 55mm, packaged 5 per zip-lock bag

MIA quad-helix arch, maxillary curved size 1-3 and maxillary straight size 1-4





Orthodontic Clasps

Ball clasps, .024 to .040, packaged 1,000 per zip-lock bag

Adams clasps, regular in 6mm to 12mm or stiff from 5mm to 12mm, packaged 100 or 1,000 per zip-lock bag

Made from high quality stainless steel



Miscellaneous

Annealed mini stop tubes, three sizes, packaged 100 per bag

AMC0165X080CRIMP

For Round up to .016"

AMC028X080CRIMP

For Rectangular up to .016" X .022",

Square up to .019" x .019" and Round up to .022"

AMC036X080CRIMP

For Rectangular up to .021 x .025",

Square up to .020" x .020"




Power pin, packaged 100 per zip-lock bag

Connected crimp stops; length of 10 stops

Power hooks, left and right, made from 316L stainless steel, packaged 10 per zip-lock bag





Straight and Cut and Spooled Lengths

Made from the same high-quality materials
as our archwires

Ideal for use in labs for creation of custom appliances

Straight lengths are packaged in tubes

Stainless steel in silver tone or gold tone, 14" long straight
and cut or 15 foot or 36 foot lengths on a bobbin, wire
sizes:

.012	.016 x .016	.018 x .022	.024	.036
.014	.016 x .022	.018 x .025	.025	.040
.016	.016 x .025	.019 x .019	.028	.045
.018	.017 x .017	.019 x .025	.030	
.020	.017 x .022	.021 x .025	.032	
	.017 x .025	.022 x .028		
	.018 x .018			

1x3 stainless steel twist in silver tone or gold tone,
14" long straight and cut, wire sizes:

.015	.0175	.0195	.018 x .025
------	-------	-------	-------------

BOBCAT 1x3 gold tone stainless steel twist,
15" long straight and cut, wire sizes:

.015	.0175	.0195	.0215
------	-------	-------	-------

Coaxial stainless steel in silver tone is available in
14" long straight and cut or 30 foot lengths on a
bobbin, wire sizes:

.0155	.0175	.0195	.0215
-------	-------	-------	-------

Nickel titanium, 14" long straight and cut or 15 foot
or 36 foot lengths on a bobbin, wire sizes:

.012	.014	.016	.018	.020
------	------	------	------	------

Beta titanium, 14" or 15" long straight and cut,
wire sizes:

.016	.017 x .025	.021 x .025
.018	.018 x .018	.027
.016 x .016	.018 x .025	.032
.016 x .022	.019 x .025	.036

Elgiloy, 14" long straight and cut, wire sizes:

.016 x .016	.016 x .022	.017 x .025	.018 x .025				
.019 x .025	.021 x .025	.028	.030	.032	.036	.040	

Auxiliaries Packaging Options

Straight and Cut Lengths:
Bulk packed in clear tube
10, 25 or 50 per clear tube


Spooled Lengths:
On bobbins

Long Ligatures:
100 or 1000 per tube

Short Ligatures:
100 or 500 per zip-lock bag

Nickel Titanium Springs:
7" lengths 3 per tube
Longer lengths on bobbin

Springs and Miscellaneous Auxiliary Items:
Packed in zip-lock bag (quantities vary)



*All orthodontic products are manufactured in
accordance with the United States Food and Drug
Administration Good Manufacturing Practices.*

Corporate Headquarters
75 Winchell Road
New Britain, CT 06052
Phone: 860-224-1349
FAX: 860-827-9982

Singapore Location
Phone: +65-6743-3374

acmemonaco.com



Linked in

©2025 Acme Monaco
All Rights Reserved
Ver. 02/2025