


BIOVis Multi-force ARCHWIRES

Acme Monaco's **BIOVis Multi-force ARCHWIRES** achieve superior clinical results by properly activating forces that are tailored to the different teeth within the dental arch. Unlike some competitive multi-force wires, Acme Monaco offers the unique ability to maintain similar clinical forces even when the material cross section changes.

Our archwire provides the optimum forces needed to torque, level, and properly align all in one solution, providing your customers the necessary tool to minimize wire changes and reduced chair time.

Acme Monaco's **BIOVis Multi-force ARCHWIRES** a high performance wire at a competitive price.




Active forces increase from anterior to posterior

Anteriors – Gentle forces in the anterior region provided by heat-active wire movement guarantees patient comfort and prevents root resorption.

Bicuspids – By providing a moderate force to this region, allows the proper biological force to be activated on the larger bicuspid.

Posteriors – Heavier force is provided in the molar region to effectively align the first and second molars.

Clinical working force at 1.5 mm
Compare to BioForce®



Product Launch
Q4 2014