

Are you seeing
a demand for a
***broader
archform?***


Xpand-form™
for a wider smile



acmemonaco.com

Xpand-form™

for a wider smile



Xpand-form is wider and with a flatter curve in the anterior section than our natural archform and it maintains expanded archform post expansion.

Xpand-form is offered in upper and lower arch forms and in all wire sizes and types with custom packaging to suit your needs.



Acme Monaco Corporation
75 Winchell Road
New Britain, CT 06052

Telephone: 860-224-1349
Email: acmecorp@acmemonaco.com
Website: acmemonaco.com